



One-stop Solution for AAV Gene Delivery



- · Online, one-stop vector design platform
- VectorBuilder has designed more than 550,000 vectors for global user
- Over 400,000 vector components avaliable
- Large variety of AAV vector expression systems


Ultra-purified Virus Packaging

- Up to 5x10¹⁴ GC per batch
- 18 AAV serotypes aviailable
- Single-stranded and self-complementary AAV
- Ready-to-ship and custom AAV empty

- directed evolution
- · One and two-step packaging approaches for high capsid diversity
- · Capsid screening and

AVAV Capsid Evolution Capsid library construction and

characterization in vivo and in vitro



AAV Blodistribution Profiling

- Mouse, rat, and non-human primate (NHP) models available
- Genome copy number analysis by multiple methods
- · Multiple routes of administration

CMP AAV Manufacturing

- Up to 10¹⁷ GC per batch
- GMP quality system for the entire process
- Analytical and process development